MANUAL ESTETOSCÓPIO BIC

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Congratulations on the purchase of your new BIC Stethoscope

GENERAL

The BIC Stethoscope is an instrument indicated for non-invasive auscultation of cardio-respiratory noises internal to the organism. The sounds are obtained by contacting the diaphragm with the region of interest. These sounds, obtained through the diaphragm, are transmitted through the sound chamber, transmitter tube, PVC tube and earphones to the operator ear. The binaural earphone allows you to hear the sound simultaneously with both ears.

CONSERVATION AND STORAGE

- Protect the stethoscope from water, moisture, impacts, falls, sunlight, heat and fire.
- Clean the stethoscope with a soft cloth, slightly moistened with water, neutral or 70% alcohol. The use of corrosive products, detergents or polishes is prohibited.
- The tube must not deform, be stretched, have sharp edges or have contact with sharp materials.
- Do not drop it and avoid shocks.
- Recommended storage temperature: - 20 ° C to 45 ° C.
- Recommended relative air humidity: maximum 80%.

HOW TO USE

Standard Adult and Pediatric Model:

- Fit the binaural earphone in your ears in order to isolate all external sounds.
- Touch the diaphragm to the area to be analyzed.
- Lightly touch the membrane a few times to see if the eartips are positioned correctly in the ears, this will produce a noise.

Adult Double Model, Pediatric Double, Rappaport and Cardiological Model:


- Fit the binaural earphone in your ears in order to isolate all external sounds.
- Lightly touch the membrane a few times to see if the eartips are positioned correctly in the ears, and which diaphragm is in the position of use, this will produce a noise.
- Touch the diaphragm to the area to be analyzed.
- Rotate the handset 180 ° to change the diaphragm to be used, for cardiac auscultation or pulmonary auscultation.

COMPOSITION AND FORMS OF PRESENTATION

Standard model


- The Standard type stethoscope has a type of headphone:
- Simple closed headphone with diaphragm.
 - 1) Angle in chromed metal and / or treated
 - 2) PVC "Y" tube
 - 3) Pair of PVC and / or silicone Olives (depending on model)

- 4) Steel spring
- 5) Headphone

Double Model

- The double model stethoscope has two types of headphones:
- Closed with rigid plastic membrane (diaphragm): used to detect low frequency sounds with greater distinction.
 - Open (bell): used to more accurately locate the sounds obtained.
 - Rotate the handset 180 ° to change the diaphragm to be used.
 - 1) Angle in chromed metal and / or treated
 - 2) PVC "Y" tube
 - 3) Pair of PVC and / or silicone Olives (depending on model)
 - 4) Steel spring
 - 5) Headphone with Isolating Ring

Rappaport model

- The Rappaport stethoscope has five different types of headphones:
- Closed type with a rigid plastic membrane (diaphragm)
 - Large: used to detect low-frequency heartbeat sounds with greater distinction.
 - Small: child use
 - Open Type: (extension bell)
 - Large: used to more accurately locate the sounds obtained at localized points: sounds with low and medium frequency.
 - Medium: Small and limited points (eg between the ribs) are concentrated to detect sounds of medium and low intensity.
 - Small: child use
 - Rotate the handset 180 ° to change the diaphragm to be used.
 - 1. Chromed metal angle
 - a. Chrome steel spring
 - b. Soft PCV Olives
 - 2. Double PVC pipe
 - c. Metal clips
 - 3. Headphone - Threaded chrome-plated brass base that allows changing the cannula or membrane.
 - d. Central rotating tube in chromed brass. It has two outputs for connection with double tube

Included:

- Adult threadable ring
- Children's threadable ring
- Pair of rigid threadable eartips
- Pair of soft threadable eartips
- Pair of threadable silicone eartips


Cardiological Model

The double model stethoscope has two types of headphones:
- Closed with rigid plastic membrane (diaphragm): used to detect low frequency sounds with greater distinction.
- Open (bell): used to more accurately locate the sounds obtained.
- Rotate the handset 180 ° to change the diaphragm to be used.
1) Stainless steel angle
2) PVC "Y" tube
3) Pair of PVC and / or silicone Olives (depending on model)
4) Steel spring
5) Matte stainless steel earpiece
a) Threaded matte stainless steel base
b) PVC and / or resin / fiber membrane
c) PVC ring

Follow:

- A pair of threadable Olives;
- A PVC and / or resin / fiber membrane

WARRANTY

Your BIC product, when new, is guaranteed against original manufacturing defects.
The warranty period starts on the date of purchase of the factory or its authorized distributors.
CBEMED's obligation is limited to the repair or replacement of components found to be defective by CBEMED within the warranty period. These warranties apply to the original purchaser and cannot be distributed or transferred to third parties. This warranty does not apply to any damage or failure of the product considered by CBEMED, in its analysis, to be due to misuse, accidents (including damage to transport), falls, neglect, improper maintenance, violation, modification or repair by persons unauthorized.
These express warranties take precedence over any and all other warranties, express or implied, including warranties of merchantability and fitness for a particular purpose, and no other person is authorized to assume on BIC's behalf any other warranties connected with the sale of the product. CBEMED is not responsible for any loss or damage, whether direct, incidental or consequential, resulting from the breach of any express warranty, except as provided herein.
The manufacturer reserves the right to replace the defective product with a new one, if deemed necessary, being the sole judgment of the manufacturer. Repairs carried out within the warranty period do not extend the warranty period.

* Standard BIC, Double BIC and Rappaport BIC stethoscope is guaranteed for 03 (three) years
* BIC Master Double Stethoscope is guaranteed for 03 (three) years
* Cardiological Stethoscope and / or belonging to the Eternity BIC Line is guaranteed for 05 (five) years

Registration with the Ministry of Health: 80540449002
Registration with the Ministry of Health: 80540449003

REGULATION

- Directive 93/42 / EEC as amended by Directive 2007/47 / EC
- Directive 93/42 / CEE amended by Directive 2007/47 / EC
- FDA - Food and Drug Administration
- ISO 13485

All warranties are valid from the date of actual delivery of the products.
For greater agility and in order to avoid improper freight, in the case of more careful analysis, CBEMED indicates that the product is sent by the customer and the manufacturing defect is confirmed, CBEMED will make a refund upon proof of value.

The reseller must complete:

Name: _____

Address: _____

Length: _____

ZIP CODE: _____

City: _____

UF: _____

Invoice No.: _____

Date ____ / ____ / ____

This warranty term must be filled out by the dealer at the time of sale. If there is a manufacturing defect, it must be presented together with the purchase invoice.

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